

**Midmark Table Bracket Part Numbers for
GS Exam Light IV, GS 300 and GS 600**

The Midmark brackets are used in conjunction with our wall/table mount.

Manual Exam Tables

Midmark part number: 9A463001

Midmark exam tables: **204, 222, 223, 604, 622, 623 and 625**

The new bracket will mount on the following exam tables:

- Ritter 204 manual exam tables
- Ritter 222/223 Barrier-Free exam tables (note: models manufactured since 2006 with removable upholstery)
- Midmark 604 manual exam tables
- Midmark 622, 623 and 625 Barrier-Free exam tables (note: models manufactured since 2006 with removable upholstery)

Power Procedure Tables

Midmark Part Number: 9A464001

* Midmark Procedure/Exam tables: **230,630, 641, 646, 647**

- Ritter 230 Power Procedure Table
- Midmark 630 Barrier-Free™ Power Procedures Table
- Midmark 641 Barrier-Free™ Power Procedures Chair
- Midmark 646 Basic Power Podiatry Procedures Chair
- Midmark 647 Barrier-Free™ Power Podiatry Procedures Chair

* Please note: These part numbers apply to the Midmark 630, 641, 646 and 647, and Ritter 230.

Most customers order procedure tables with the base rails, where the new light is mounted.

If they are ordering a **new** procedure table, they can order the base rails factory installed:

- **9A205003** Base rails, plant installed

If they are ordering base rails to **retrofit** on site, they would order:

- **9A205004** Base rails, field installed

Ordering Information for both brackets and Accessory Base Rails:
Distribution or directly through Midmark @ 1-800-MIDMARK