



Welch Allyn MacroView™ Veterinary Otoscope

FOR VETERINARIANS



A Difference You've Got to See to Believe

Get a wider, deeper, clearer view of the ear with the Welch Allyn MacroView™ Veterinary Otoscope, the latest advance to the device that helps you perform one of your most common procedures. Using advanced technology in an innovative design, the MacroView™ Veterinary Otoscope gives you a large, sharp, nearly full view of the tympanic membrane—even in large dogs—virtually eliminating the need for panning. And with an easy-to-focus eyepiece, you can even adjust focus for different animals and ear canal lengths, meaning you get crystal-clear views of the ear canal at any point—allowing you to more effectively diagnose and treat your patients faster than ever before.

The Welch Allyn MacroView™ Veterinary Otoscope gives you:

- More detail with better clarity—approximately 27% greater magnification than with traditional otoscopes
- Approximately 35% greater field of view than with traditional otoscopes
- The ability to adjust focus for variable ear canal lengths and can correct for farsighted eyes
- Tip Grip for secure fastening of ear specula, along with ejection feature for easy removal
- Built-in throat illuminator

WelchAllyn®

Advancing Frontline Care™


The Welch Allyn MacroView™ Veterinary Otoscope Gives You a Bigger, Better, Sharper View

FOR VETERINARIANS




The new instrumentation tip gives you the advantage of both magnification and focus for easier instrumentation.

INCREASED FIELD OF VIEW




**MacroView™
Veterinary Otoscope**

Get a nearly full, crystal-clear view of the tympanic membrane—even in large canines!


“The wider angle of view and ability to focus the MacroView™ Veterinary Otoscope means better evaluation by my doctors and better outcomes for our patients. Finally, I’m able to fully visualize the ear canal and tympanum of my pet patients!”

—Dr. Ernest Ward, Jr.
Seaside Animal Care, Calabash, N.C.



Traditional Otoscope

LONGER, ADJUSTABLE FOCAL LENGTH



Ordering Information

23862 MacroView Veterinary Otoscope with Specula

Accessories

- 23834 4 mm MacroView Veterinary Reusable Specula
- 23835 5 mm MacroView Veterinary Reusable Specula
- 23837 7 mm MacroView Veterinary Reusable Specula
- 23842 MacroView Veterinary Reusable Specula Set
- 23847 MacroView Veterinary 7 mm Instrumentation Specula
- 06500-U Replacement Lamp

Diagnostic Sets featuring the MacroView Veterinary Otoscope are also available.

To find out how the new Welch Allyn MacroView™ Veterinary Otoscope can give you bigger, sharper views of the ear, go to www.welchallyn.com/dvm or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA
(p) 1.800.535.6663 (f) 1.315.685.3361 www.welchallyn.com/dvm

WelchAllyn®

Advancing Frontline Care™